



# LIPOSAT®, VIBRASAT® AND VACUSAT® POWER

## SYSTEMS FOR AESTHETIC MEDICINE

The technically leading system, a anaesthetic infiltration pump (Liposat® power) and a suction system (Vacusat® power) with vibrating handpiece (Vibrasat® power) from Möller Medical, is widely used by beauty clinics and medical practices. The fully autoclavable hand-held device and the cannula for fast, powerful axial movements facilitate the doctor's work.

- Gentle, homogeneous liposuction
- Less effort needed, allowing the doctor to work with greater precision
- Low anaesthetic risk as patient stays conscious
- Low loss of blood
- Minimal pain after the (outpatient) operation
- CE-certified



### Liposuction

Liposuction is the most common surgical procedure in plastic surgery. There are estimated to be as many as 400,000 such operations each year worldwide, tendency rising. Tumescence local anaesthesia with vibration liposuction is a gentle form of treatment. This involves injecting between two and three litres of salt solution containing an anaesthetic into the fatty tissue. Once it has had time to take effect, a cannula connected to a liposuction device is introduced through small cuts in the skin. The doctor moves the cannula vigorously through the tissue – or uses an electric vibration drive unit – to remove the fat.



## LIPOSAT<sup>®</sup>, VIBRASAT<sup>®</sup> AND VACUSAT<sup>®</sup> POWER

### SYSTEMS FOR AESTHETIC MEDICINE

#### Vibrasat<sup>®</sup> power

Vibrasat<sup>®</sup> power is the technologically leading vibrating device for liposuction.

- Up to 5000 strokes of 2.8 mm per minute
- Very powerful; about four times more powerful than its predecessor product Vibrasat<sup>®</sup>
- Fully autoclavable thanks to the use of special materials for the drive motor and gears
- Anatomically shaped handpiece with just the right centre of gravity and low weight for ease of working, yet sufficiently heavy to cushion the hand adequately against the vibrations
- Medically approved device bearing the CE mark
- Very quiet-running
- Quality product made in Germany

#### Liposat<sup>®</sup> power / Vacusat<sup>®</sup> power

Equipped with the Liposat<sup>®</sup> power infiltration pump, the Vacusat<sup>®</sup> power suction unit and the Vibrasat<sup>®</sup> power vibrating handpiece, plastic/aesthetic surgeons and dermatologists can perform liposuction treatments quickly and easily at their practices.

- Preselection options on the infiltration pump save time: instead of having to inject the saline solution gradually, the doctor can specify the volume and delivery rate of the Liposat<sup>®</sup> power and then concentrate on other preparatory activities.
- Heating dishes keep two bags of preheated tumescence solution at a constant temperature. Safe and comfortable for the patient.
- Pressure up to 3 bar; delivery rate up to 300 ml/min
- Liposat<sup>®</sup> power, the consumables and accessories are medically licensed for infiltration and bear the CE mark.
- Quiet-running
- Quality product made in Germany



A Product of

Möller Medical GmbH & Co. KG  
Wasserkuppenstrasse 29 - 31  
D-36043 Fulda

Tel. +49 (0) 661 941 95 0

Fax +49 (0) 661 941 95 90

E-mail: [info@moeller-medical.com](mailto:info@moeller-medical.com)

[www.moeller-medical.com](http://www.moeller-medical.com)